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

INSTRUCTIONS FOR USE

for optimed percutaneous introducers / percutaneous introducer sets

Please Note

The product is sterile packaged (EO sterilisation), designed for single use and must not be re-sterilized. The product must be stored in a cool, dry place and protected from light. This product must not be used if either product or packaging are damaged. **Only a doctor familiar with its application may use this product. This medical product is designed for single use only. It must not be reprocessed, cleaned, disinfected or sterilized!** After reprocessing, the product is no longer suitable for use and may constitute a biological risk for both patient and user: owing to the surface composition and design of the product, full effectiveness of cleaning after use cannot be guaranteed and contamination cannot be ruled out if the product is re-used. Re-sterilisation may impair the surface composition and material properties of the plastics used. In this case the original function of the product cannot be guaranteed.

Description / Structure

Percutaneous Introducers	Percutaneous Introducer Sets
① Sheath	① Sheath
② Dilator	② Dilator
③ Removable 3-way stopcock	③ Removable 3-way stopcock
④ (Removable) hemostatic valve	④ (Removable) hemostatic valve
	⑤ Guide wire
	⑥ Angiography needle

Phthalates

This medical product may contain phthalates. Please check the label on the packaging for information. Phthalates (plasticizers) are used to increase the flexibility of the material. The product has been designed and produced in such a way as to minimize risks arising from substances released by the product. This also applies to the use of the product with children and pregnant or breastfeeding women.

Note

Indications, contraindications and application of the above product must always comply with current medical standards.

Indications

Facilitation of introduction of a variety of diagnostic and therapeutic instruments into the peripheral vascular system

Absolute Contra-Indications

- No safe percutaneous access possible
- Non-treatable coagulation disorders

Possible Complications

- Only doctors familiar with the possible complications may apply this product. Complications may occur at any time during or following the procedure. Possible complications include:
 - Dislocation
 - Tissue damage
 - Inkustation / occlusion
 - Hemorrhage/hematoma
 - Vascular perforation / Dissection
 - Iatrogenic infections / sepsis
 - Thromboses / embolism
 - Pain
 - Death
- The percutaneous introducer/introducer set is not suitable for use in the central circulatory system (as per EC regulation 93/42/EEC, annex IX, points 1.7. and 1.8.).

Preparation

- Check packaging for damage. If the packaging is damaged, the product must not be used.
- Remove product from the packaging.
- Check product for damage. A damaged product must not longer be used.
- Ensure that products not belonging to the set are compatible.

Registration

Please see the **optimed** brochures and / or product labels for precise details of all products listed here. The package includes a removable label where product details are registered

Notes on Correct Use

- Only introduce the percutaneous introducer under fluoroscopic control.
- Do not introduce the sheath without the guide wire and inner dilator.
- Do not cut the guide wires, introducer sheath and dilators. Changes to the set may impair its function.
- Care must be taken not to damage the sheath when penetrating, suturing and cutting tissue in the vicinity of the puncture point.
- Before repeating the introduction of a medical instrument, aspirate the contents of the sheath through the side arm of the valve and flush in line with medical standards.

Use

Take care to keep the sheath aseptic in compliance with current medical standards.

Preparing the Patient

- Prepare puncture site in compliance with current medical standards.
- Locally anaesthetise the puncture site in accordance with current medical standards.

Preparing the Sheath

- Flush the dilator, introducer sheath and side arm with heparinized physiological saline solution.
- Introduce the dilator into the sheath through the hemostatic valve and lock.

Puncture

- Introduce the puncture needle into the respective vessel and set the guide wire.
- When the guide wire is correctly positioned, remove the puncture needle (Note: hold the guide wire firmly).
- Now use the guide wire to insert the sheath and dilator into the lumen of the vessel by slightly rotating the device.
- Check the position of the sheath by withdrawing the guide wire and dilator until blood flows through the side arm of the flushing unit.
- Remove the dilator using the guide wire
- Fix the sheath in position
- Replace the short guide wire with a standard guide wire
- Now the desired instruments can be positioned using the hemostatic valve.

Note

- If resistance is encountered do not extend or withdraw the sheath. Likewise, if resistance is encountered do not extend or withdraw any interventional instruments (e.g. guide wire, catheter) through the dilator and/or sheath.
- Continuation of intervention may seriously injure vessels and/or damage the components of the interventional instruments.
- Before proceeding, it is essential to remove the cause of resistance (using fluoroscopic control if necessary).

Removing the Sheath

- Remove the catheter from the sheath
- Slowly remove the sheath in parallel position to the skin
- Dress puncture wound according to current medical standards

Further procedure according to current medical standards.

Recommended Medication

Preoperative, intraoperative and postoperative medication should be initiated in compliance with the current medical standards.

Medical products may constitute a biological risk after use. They should only be used and disposed of in accordance with statutory regulations and acknowledged practice.

Limited Guarantee

optimed guarantees that its products are manufactured with the greatest possible care. THIS IS THE ONLY VALID GUARANTEE AND SUBSTITUTES ALL OTHER GUARANTEE DECLARATIONS WHICH HAVE BEEN GIVEN. It should be noted that no product is always effective under all circumstances; this is due to the biological differences of the persons to be treated. Components of **optimed** sets as well as individual **optimed** products are interchangeable. Because using individual **optimed** products or sets with non-**optimed** products, the user must ensure the compatibility of the individual products with the specific application. **optimed** has no influence on how the product is used, the diagnosis of the patient, nor how the product is handled outside of the company. **optimed** cannot guarantee the effectiveness nor the absence of complications. Therefore, **optimed** assumes no liability for damages and costs. **optimed** will replace products which have a defect for which **optimed** is responsible. **optimed** cannot be held liable for consequential damage of any kind which was caused by resterilising or reprocessing the product. **optimed** employees are not entitled to change the above terms and conditions, to extend the liability, or to assume additional obligations concerning the product. Products are subject to alterations.

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